

## PREASSURE ACTIVATED SAFETY LANCETS

### PRODUCT INFORMATION AND DIRECTIONS FOR USE

#### Features and Benefits

Safety Lancets are colour-coded per gauge for easy recognition | There are different penetration depths and needle gauges available for different skin types | Easy to use just twist and click | Triple sharpened needle tips, increasing patient comfort | Ergonomic design for comfortable grip | Automatic retraction of needle or blade to eliminate the possibility of needle stick injuries.

#### Directions for Use:

Clean test site prior to use | While firmly holding the body of the device, twist off the protective cap | Place the device on the test site and push firmly against the skin | After use dispose the safety lancet in the appropriate sharp's container.



\*\*All pictures shown are for illustration purpose only. Actual colours may vary.\*\*

**CE** 0197  
Medical Device Class IIa



ISO 13485

**STERILE** R

STERILIZED BY GAMMA RADIATION